



Recombinant Arabidopsis thaliana SPX domain-containing protein 4 (SPX4)

Product Code	CSB-MP850081DOA
Abbreviation	SPX4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94A21
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MKFGKEFRTH LEETLPEWRD KFLCYKPLKK LLKYYPYISA DFGPANSDHN DSRPVFADTT NISSAADDGG VVPGVRPSED LQGSFVRILN DELEKFNDFY VDKEEDFVIR LQELKERIEQ VKEKNGEFAS ESEFSEEMMD IRRDLVTIHG EMVLLKNYSS LNFAGLVKIL KKYDKRTGGL LRLPFTQLVL HQPFTEPL TRLVRECEAN LELLFPSEAE VVESSAVQA HSSSHQHNSP RISAETSSTL GNENLDIYKS TLAAMRAIRG LQKASSTYNP LSFSSLLQNE DDETVTAENS PNSGNKDDSE KEDTGPSH
Source	Mammalian cell
Target Names	SPX4
Protein Names	Recommended name: SPX domain-containing protein 4 Alternative name(s): Protein SPX DOMAIN GENE 4 Short name= AtSPX4
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.