



Recombinant *Oryza sativa* subsp. *japonica* Tubby-like F-box protein 2 (TULP2)

Product Code	CSB-MP850095OFG
Abbreviation	TULP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94DT9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MVPWRRSSSS SSAPSSRPAR RPARTNARVS PDVSSSELSPL AGEEGAGEER WSALVPDLLA DILRCVEAGS ERWPPRRDVV ACASVCRRWR DVAVAVVQPP LESGKITFPS SLKQPGPRDA PMQCFIKRKN KNSTFFLYLG LTQELTDDEK FLLAARRCRR GLHKEYAITI NSDGLFHGSQ SCVGNLKS NF TGTKFTIRDW QPPYEGAKAF SSRSGRWFGN KHRCLVSTG DVEVGEVSYK YSLLRPRGPR RMSCSVQCPV LKGTAVDPQD GKRLSNSIPS SLVLNSKVPS WHEHLQCWCL NFHGRVMVAS VKNFQLIAPV EPGEPSDKTV VLQFGKIDDD VFTMDYRQPL SAFQAFAICL SNFGTKLA
Source	Mammalian cell
Target Names	TULP2
Protein Names	Recommended name: Tubby-like F-box protein 2 Short name= OsTLP2 Alternative name(s): Tubby-like F-box protein 12 Short name= OsTLP12
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.