



Recombinant *Oryza sativa* subsp. japonica Eukaryotic translation initiation factor 3 subunit K (TIF3K1)

Product Code	CSB-BP850101OFG
Abbreviation	TIF3K1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q94HF1
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASEQAAESY TVEELVAVNP YNPDILNDLE GFVNDQVSNQ TYNLDANLSL LRLYQFEPER LSVQIVSRIL IKALMAMPGP DFSLCLFLIP EHVQMEEQFK TLIVLSHYLE TARFRQFWDE ASKNRNILDV VPGFEQAIQS YAIHVLSLTY QKVPRPVLAE AINIEGLALD KFLEHHIANS GWVIEKGARS QLIVLPRNEF NHPELKKNTA ETVPFHEHVTR IFPVLS
Source	Baculovirus
Target Names	TIF3K1
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit K Short name= eIF3k Alternative name(s): eIF-3 p25
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.