




Recombinant *Macaca fascicularis* T-cell surface glycoprotein CD3 gamma chain (CD3G), partial

Product Code	CSB-YP850179MOV
Relevance	The CD3 complex mediates signal transduction.
Abbreviation	Recombinant Cynomolgus monkey CD3G protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95LI7
Alias	T-cell receptor T3 gamma chain CD_antigen: CD3g
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (Cynomolgus monkey)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	QSFEENRKLNVYNQEDGSVLLTCHVKNTNITWFKEGKMGIDILTAHKNKWNLGS NTKDPRGVYQCKGSKDKSKTLQVYYRMCQNCIELNAAT
Research Area	Immunology
Source	Yeast
Target Names	CD3G
Protein Names	Recommended name: T-cell surface glycoprotein CD3 gamma chain Alternative name(s): T-cell receptor T3 gamma chain CD_antigen= CD3g
Expression Region	23-113aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	12.5kDa
Protein Length	Extracellular Domain

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.