



# Recombinant Cat Interleukin-18 (IL18)

<b>Product Code</b>	CSB-MP850184CA
<b>Abbreviation</b>	IL18
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q95M33
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Felis catus (Cat) (Felis silvestris catus)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	YFGKL EHKLSILRNL NDQVLFINQG DQPVFEDMPD SDCTDNAPRT EFIIYMYKDS LTRGLAVTIS VNYKTMSTLS CENKIISFKE MSPPESINDE GNDIIFQRS VPGHDDKIQF ESSLYKGYFL ACEKEKDLFK LILKKKDENG DKSIMFTVQN KN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IL18
<b>Protein Names</b>	Recommended name: Interleukin-18 Short name= IL-18 Alternative name(s): Interferon gamma-inducing factor Short name= IFN-gamma-inducing factor Interleukin-1 gamma Short name= IL-1 gamma
<b>Expression Region</b>	36-192
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a proinflammatory cytokine. This cytokine can induce the IFN-gamma production of T cells. The combination of this cytokine and IL12 has been shown to inhibit IL4 dependent IgE and IgG1 production, and enhance IgG2a production of B cells. IL-18 binding protein (IL18BP) can specifically interact with this cytokine, and thus negatively regulate its biological activity.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.