



Recombinant Gorilla gorilla gorilla Baculoviral IAP repeat-containing protein 8 (BIRC8)

Product Code	CSB-EP850185GGZ-B
Abbreviation	BIRC8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95M71
Product Type	Recombinant Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Purity	>85% (SDS-PAGE)
Sequence	MTGYEAWLIT FGTWMYSV NK EQLARAGFYA IGQEDKIQCF HCGGGLANWK PKEDPWEQHA KWYPGCKYLL EEKGHEYINN IHLTRSLEGA LVQTTKKTPS LTKRISDTIF PNPMLQEAIR MGFDFKDIKK IMEEKIQTSG SNYKTLEVLV ADLVSAQKDT TENESNQTSL QREISPEEPL RRLQEEKLCK ICMDRHIADV FIPCGHLVTC KQCAEAVDRC PMCNAVIDFK QRVFMS
Source	E.coli
Target Names	BIRC8
Protein Names	Recommended name: Baculoviral IAP repeat-containing protein 8 Alternative name(s): Inhibitor of apoptosis-like protein 2 Short name= IAP-like protein 2 Short name= ILP-2
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.