



Recombinant PCI domain-containing protein 2 homolog (C27F2.10)

Product Code	CSB-MP850197CXY
Abbreviation	C27F2.10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q95QU0
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPSVRSIEDY FGQIESLLYR QDWTNGEKIS KFVSTYDEHA QEPFMHIEAY GSRSKRCRVS EDEVFDEIVC LHLHVLYNIH VAQDLITAQS TQIQIIQLFN KEILQKRKDE NWFLPIFYRL CTDLRWLSKG AEACVSGDDE GDSNANSFFE SAAKAITECY RTCVSDVHAE EGTTKKVAML NMTNQLFQIY FQINKLNLLK PLIRAIDNCG SLYHDFLMSD KVAINYFLGR KAMFDADLNL AEKSLLYAFR NCPADSMSNK RKILIYLIPV KMFLGHMPTS QLLHEYRLDE FQDVVAGVKD GNLAQLDGAL AANEAFFIKC GIFLMLEKLR MITFRTLFFK VSQIVGTAQI PLDAFQTALR FVGVTDVDMD ELECIANLI ASKKIKGYLS HQHQKLVISK MNAFPTLSGV SSN
Source	Mammalian cell
Target Names	C27F2.10
Protein Names	Recommended name: PCI domain-containing protein 2 homolog Alternative name(s): CSN12-like protein
Expression Region	1-413
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.