



# Recombinant Human Interleukin-22 receptor subunit alpha-2 (IL22RA2)

<b>Product Code</b>	CSB-BP850247HU
<b>Abbreviation</b>	IL22RA2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q969J5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TQSTHESLK PQRVQFQSRN FHNILQWQPG RALTGNSSVY FVQYKIMFSC SMKSSHQKPS GCWQHISCNF PGCRTLAKYG QRQWKNKEDC WGTQELSCDL TSETSDIQEP YYGRVRAASA GSYSEWSMTP RFTPWWETKI DPPVMNITQV NGSLLVILHA PNLPHYRYQKE KNVSIEDYYE LLYRVFIINN SLEKEQKVYE GAHRAVEIEA LTPHSSYCVV AEIYQPMLDR RSQRSEERCV EIP
<b>Source</b>	Baculovirus
<b>Target Names</b>	IL22RA2
<b>Protein Names</b>	Recommended name: Interleukin-22 receptor subunit alpha-2 Short name= IL-22 receptor subunit alpha-2 Short name= IL-22R-alpha-2 Short name= IL-22RA2 Alternative name(s): Cytokine receptor class-II member 10 Cytokine recept
<b>Expression Region</b>	22-263
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This protein is a soluble class II cytokine receptor. This protein has been shown to specifically bind to interleukin 22 (IL22), block the interaction of IL22 with its cell surface receptor, and thus inhibit IL22 activity. This protein functions as an IL22 antagonist, and may be important in the regulation of inflammatory response. Three alternatively spliced transcript variants encoding distinct isoforms have been described.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.