



Recombinant Human Peptidylprolyl isomerase domain and WD repeat-containing protein 1 (PPWD1)

Product Code	CSB-YP850265HU
Abbreviation	PPWD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BP3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AAESGSDFQ QRRRRRRDPE EPEKTELSER ELAVAVAVSQ ENDEENEERW VGPLPVEATL AKKRKVFLEFE RVYLDNLPSA SMYERSYMHR DVITHVVCTK TDFIITASHD GHVKFWKKIE EGIEFVKHFR SHLGVIESIA VSSEGALFCS VGDDKAMKVF DVVNFDMINM LKLGYPFGQC EWIIYCPGDAI SSVAASEKST GKIFIYDGRG DNQPLHIFDK LHTSPLTQIR LNPVYKAVVS SDKSGMIEYW TGPPHEYKFP KNVNWEYKTD TDLYEFAKCK AYPTSVCFSP DGKKIATIGS DRKVRIFRFV TGKLMRVFDE SLSMFTLQQ MRQQLPDMEF GRRMAVEREL EKVDAVRLIN IVFDETGHFV LYGTMLGIKV INVETNRCVR ILGKQENIRV MQLALFQGIA KKHRAATTIE MKASENPVLQ NIQADPTIVC TSFKKNRFYM FTKREPEDTK SADSDRDVFN EKPSKEEVMA ATQAEGPKRV SDSAIHTSM GDIHTKLFVP ECPKTVENFC VHSRNGYYNG HTFHRIIKGF MIQTGDPTGT GMGGESIWGG EFEDEFHSTL RHDRPYTL SM ANAGSNTNGS QFFITVVPTP WLDNKHTVFG RVTKGMEVVQ RISNVKVNPK TDKPYEDVSI INITVK
Source	Yeast
Target Names	PPWD1
Protein Names	Recommended name: Peptidylprolyl isomerase domain and WD repeat-containing protein 1 EC= 5.2.1.8 Alternative name(s): Spliceosome-associated cyclophilin
Expression Region	2-646
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.