



Recombinant Human MOB kinase activator 3A (MOB3A)

Product Code	CSB-MP850266HU
Abbreviation	MOB3A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96BX8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSNPFLKQVF NKDKTFRPKR KFEPGTQRFE LHKKAQASLN AGLDLRLAVQ LPPGEDLNDW VAVHVVDFFN RVNLIYGTIS DGCTEQSCPV MSGGPKYEYR WQDEHKFRKP TALSAPRYMD LLMDWIEAQI NNEDLFPTNV GTPFPKNFLQ TVRKILSRLF RVFVHVYIHH FDRIAQMGE AHVNTCYKHF YYFVKEFGLI DTKELEPLKE MTARMCH
Source	Mammalian cell
Target Names	MOB3A
Protein Names	Recommended name: MOB kinase activator 3A Alternative name(s): MOB-LAK Mob1 homolog 2A Mps one binder kinase activator-like 2A
Expression Region	1-217
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.