



Recombinant Human Leucine-rich repeat-containing protein 45 (LRRC45)

Product Code	CSB-BP850272HU
Abbreviation	LRRC45
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96CN5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEEFRSYSR LCRESGAEPQ EAVLQQLHQL PRGRDLATQ SLTVETCRAL GKLLPRETLC TELVLSDCML SEEGATLLLR GLCANTVLRV LDLKGNNLRA AGAEALGKLL QQNKSIQSLT LEWNSLGTWD DAFATFCGGL AANGALQRLD LRNNQISHKG AEELALALKG NTTLQQLDLR WNNVGLLGGR ALMNCLPSNR TLWRLDLAGN NIPGDVLRV EQAMGHSQDR LTTFQENQAR THVLSKEVQH LREEKSKQFL DLMETIDKQR EEMAKSSRAS AARVGQLQEA LNERHSIINA LKAKLQMTEA ALALSEQKAQ DLGELLATAE QEQLSLSQRQ AKELKLEQQE AAERESKLLR DLSAANEKNL LLQNQVDELE RKFRCQQEQL FQTRQEMTSM SAELKMRAIQ AEERLDMEKR RCRQSLEDSE SLRIKEVEHM TRHLEESEKA MQERVQRLEA ARLSLEEELS RVKAAALSER GQAEELIKA KSQARLEEQQ RLAHLEDKLR LLAQARDEAQ GACLQQKQVV AEAQTRVSQL GLQVEGLRRR LEELQQELSL KDQERVAEVS RVRVELQEQN GRLQAEELAAQ EALREKAAAL ERQLKVMASD HREALDRES ENASLREKLR LREAELIARIR DEEAQRASFL QNAVLAYVQA SPVRTLSPPK
Source	Baculovirus
Target Names	LRRC45
Protein Names	Recommended name: Leucine-rich repeat-containing protein 45
Expression Region	1-670
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.