



Recombinant Human Protein SAAL1 (SAAL1)

Product Code	CSB-EP850288HU-B
Abbreviation	SAAL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96ER3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDRNPSPPPP GRDKEEEEEV AGGDCIGSTV YSKHWLFGVL SGLIQIVSPE NTKSSSDDEE QLTELDEEME NEICRVWDMS MDEDVALFLQ EFNAPDIFMG VLAASKCPRL REICVGILGN MACFQEICVS ISSDKNLGQV LLHCLYDSDP PTLLETSRLL LTCLSQA EVA SVWVERIQEH PAIYDSICFI MSSSTNVDLL VKVGEVVDKL FDLDEKLMLE WVRNGAAQPL DQPQEESEEQ PVFRLVPCIL EAAKQVRSEN PEWLDVYMHI LQLLTTVDDG IQAIVHCPDT GKDIWNLLFD LVCHEFCQSD DPPIILQE QK TVLASVFSVL SAIYASQTEQ EYLKIEKVDL PLIDSLIRVL QNMEQCQKKP ENSAESNTEE TKRTDLTQDD FHLKILKDIL CEFLSNIFQA LTKETVAQGV KEGQLSKQKC SSAFQNLLPF YSPVVEDFIK ILREVDKALA DDLEKNFPSL KVQT
Source	E.coli
Target Names	SAAL1
Protein Names	Recommended name: Protein SAAL1
Expression Region	1-474
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.