



# Recombinant Human Cell division cycle-associated 7-like protein (CDCA7L)

<b>Product Code</b>	CSB-EP850301HU
<b>Abbreviation</b>	CDCA7L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96GN5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MELATRYQIP KEVADIFNAP SDDEEFVGFR DDVPMETLSS EESCDSFDSL ESGKQQDVRF HSKYFTEELR RIFIEDTDSE TEDFAGFTQS DLNGKTNPEV MVVESDLSDG GKASLVSEEE EDEEEDKATP RRSRSTRSSI GLRVAQFQPT KKLANKPKDN SSSEQLFSSA RLQNEKKTIL ERKKDCRQVI QREDSTSESE DDSRDESEQES SDALLKRTMN IKENKAMLAQ LLAELNSMPD FFPVRTPTSA SRKKTVRRAF SEGQITRRMN PTRSARPPEK FALENFTVSA AKFAEEFYSF RRRKTIGGKC REYRRRHRIS SFRPVEDITE EDLENAVITV RDKIYDKVLG NTCHQCRQKT IDTKTVCRNQ GCCGVRGQFC GPCLRNRYGE DVRSAALLDPD WVCPPCRGIC NCSYCRKRDG RCATGILIH L AKFYGYDNVK EYLESLQKEL VEDN
<b>Source</b>	E.coli
<b>Target Names</b>	CDCA7L
<b>Protein Names</b>	Recommended name: Cell division cycle-associated 7-like protein Alternative name(s): Protein JPO2 Transcription factor RAM2
<b>Expression Region</b>	1-454
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.