



Recombinant Human Cell division cycle-associated 7-like protein (CDCA7L)

Product Code	CSB-MP850301HU
Abbreviation	CDCA7L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96GN5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MELATRYQIP KEVADIFNAP SDDEEFVGFR DDVPMETLSS EESCDSFDSL ESGKQQDVRF HSKYFTEELR RIFIEDTDSE TEDFAGFTQS DLNGKTNPEV MVVESDLSDG GKASLVSEEE EDEEEDKATP RRSRSTRSSI GLRVAFQFPT KKLANKPKDN SSSEQLFSSA RLQNEKKTIL ERKKDCRQVI QREDSTSESE DDSRDESEQES SDALLKRTMN IKENKAMLAQ LLAELNSMPD FFPVRTPTSA SRKKTVRRAF SEGQITRRMN PTRSARPPEK FALENFTVSA AKFAEEFYSF RRRKTIGGKC REYRRRHRIS SFRPVEDITE EDLENAVITV RDKIYDKVLG NTCHQCRQKT IDTKTVCRNQ GCCGVRGQFC GPCLRNRYGE DVRSAALLDPD WVCPPCRGIC NCSYCRKRDG RCATGILIHLL AKFYGYDENVK EYLESLQKEL VEDN
Source	Mammalian cell
Target Names	CDCA7L
Protein Names	Recommended name: Cell division cycle-associated 7-like protein Alternative name(s): Protein JPO2 Transcription factor RAM2
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.