



Recombinant Human SH2 domain-containing adapter protein D (SHD)

Product Code	CSB-EP850314HU-B
Abbreviation	SHD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96IW2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAKWLRDYLS FGGRPPPPQP PTPDYTESDI LRAYRAQKNL DFEDPYEDAE SRLEPDPAGP GDSKNPGDAK YGSPKHRLIK VEAADMARAK ALLGGPGEEL EADTEYLDPF DAQPHPAPPD DGYMEPYDAQ WVMSELPGRG VQLYDTPYEE QDPETADGPP SGQKPRQSRM PQEDERPADE YDQPWEWKKD HISRAFAVQF DSPEWERTPG SAKELRRPPP RSPQPAERVD PALPLEKQPW FHGPLNRADA ESLLSLCKEG SYLVRLSETN PQDCSLSLRS SQGFLHLKFA RTRENQVVLG QHSGPFPSVP ELVLHYSSRP LPVQGAEHLA LLYPVVTQTP
Source	E.coli
Target Names	SHD
Protein Names	Recommended name: SH2 domain-containing adapter protein D
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.