



Recombinant Human Ubiquitin-conjugating enzyme E2 E2 (UBE2E2)

Product Code	CSB-EP850386HU-B
Abbreviation	UBE2E2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96LR5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	STEAQRVDD SPSTSGGSSD GDQRESVQQE PEREQVQPKK KEGKISSKTA AKLSTSAKRI QKELAEITLD PPPNCSAGPK GDNIYEWRSST ILGPPGSVYE GGVFFLDITF SPDYPFKPPK VTFRTRIIYHC NINSQGVICL DILKDNWSPA LTISKVLLSI CSLLTDCNPA DPLVGSIIATQ YMTNRAEHDR MARQWTKRYA T
Source	E.coli
Target Names	UBE2E2
Protein Names	Recommended name: Ubiquitin-conjugating enzyme E2 E2 EC= 6.3.2.19 Alternative name(s): Ubch8 Ubiquitin carrier protein E2 Ubiquitin-protein ligase E2
Expression Region	2-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.