



# Recombinant Human Uncharacterized protein C15orf33 (C15orf33)

<b>Product Code</b>	CSB-EP850390HU-B
<b>Abbreviation</b>	C15orf33
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q96M60
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAGQRTCQRR SSRAGPGKMQ EPPKSIEEFL KFQNWWDYWPR EIHFRDDDKW SCTLKKIKED SSFVSIYTHL WENVPRIFEA LLIMESKLKE YSLILQNHTS EIFKWKSMIS ETSSYRKLER YGEFLKKYHK KKKIMLSDEM ETEKNIEGCS FTGFKANELT QLPRHLDAEQ IYLFILKAHN FDERVFKIWK THFLSEASIA LLHDSFWWWF LHKFRPDREN QDCLFDRISE SYVTLFMSIP LSRKDAFFQI YPDCLAQAIY ATFHEAFPES SYLFNDEFKE DLGNNIFLWC SGLKPQKGFW IHWKLKELST TTIHGSKKAP AKSVKERIAD SQEHISTSID FNIKILNNP RAYTLPISKE ESRLSRLATK SHYSSTGPEF NRVLFNFGGQ SPLILYYLKM HELAGISKAP KKTIKIKLTKI FQEPLPAPTY RDVIKEAKRQ FARNQKDFRI LQAKATKKPH EVKQDFEKFL HKLRSEAEIE RECVASLSSS SSSSPSTDN YNFEEEEY
<b>Source</b>	E.coli
<b>Target Names</b>	FAM227B
<b>Protein Names</b>	Recommended name: Uncharacterized protein C15orf33
<b>Expression Region</b>	1-508
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.