



Recombinant Human Uncharacterized protein C15orf33 (C15orf33)

Product Code	CSB-MP850390HU
Abbreviation	C15orf33
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96M60
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAGQRTCQRR SSRAGPGKMQ EPPKSIEEFL KFQNWDYWPR EIHFRDDDKW SCTLKKIKED SSFVSIYTHL WENVPRIFEA LLIMESKLKE YSLILQNHTS EIFKWKSMIS ETSSYRKLER YGEFLKKYHK KKKIMLSDEM ETEKNIEGCS FTGFKANELT QLPRHLDAEQ IYLFILKAHN FDERVFKIWK THFLSEASIA LLHDSFWWWF LHKFRPDREN QDCLFDRISE SYVTLFMSIP LSRKDAFFQI YPDCLAQAIY ATFHEAFPES SYLFNDEFKE DLGNNIFLWC SGLKPQKGFW IHWKLKELST TTIHGSKKAP AKSVKERIAD SQEHISTSID FNIKILNNP RAYTLPISKE ESRLSRLATK SHYSSTGPEF NRVLFNFGGQ SPLILYYLKM HELAGISKAP KKTIKIKLTKI FQEPLPAPTY RDVIKEAKRQ FARNQKDFRI LQAKATKKPH EVKQDFEKFL HKLRSEAEIE RECVASLSSS SSSSPSTDN YNFEEEEY
Source	Mammalian cell
Target Names	FAM227B
Protein Names	Recommended name: Uncharacterized protein C15orf33
Expression Region	1-508
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.