



Recombinant Human BTB/POZ domain-containing protein KCTD7 (KCTD7)

Product Code	CSB-BP850395HU
Abbreviation	KCTD7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96MP8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVVVTGREPD SRRQDGAMSS SDAEDDFLEP ATPTATQAGH ALPLLPQEFP EVVPLNIGGA HFTTRLSTLR CYEDTMLAAM FSGRHYIPTD SEGRYFIDRD GTHFGDVLNF LRSGDLPPRE RVRVAVYKEAQ YYAIGPLLEQ LENMQPLKGE KVRQAFLGLM PYYKDHLEI VEIARLRAVQ RKARFAKLKV CVFKEEMPIT PYECPLLNSL RFERSESDGQ LFEHHCEVDV SFGPWAEVAD VYDLLHCLVT DLAQGLTVD HQCIGVCDKH LVNHYYCKRP IYEFKITWW
Source	Baculovirus
Target Names	KCTD7
Protein Names	Recommended name: BTB/POZ domain-containing protein KCTD7
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.