



Recombinant Human Rab-interacting lysosomal protein (RILP)

Product Code	CSB-BP850399HU
Abbreviation	RILP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NA2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEPRRAAPGV PGWGSREAAG SASAAELVYH LAGALGTELQ DLARRFGPEA AAGLVPLVVR AELLEQA AV GPAPDSLQVS AQP AEQELRR LREENERLRR ELRAGPQEER ALLRQLKEVT DRQRDELRAH NRDLRQRGQE TEALQEQLQR LLLVNAELRH KLAAMQTQLR AAQDRERERQ QPGEAATPQA KERARGQAGR PGHQHGQEPE WATAGAGAPG NPEDPAEAAQ QLGRPSEAGQ CRFSREEFEQ ILQERNELKA KVFLIKEELA YFQRELLTDH RVPGLLLEAM KVAVRKQRKK IKAKMLGTPE EAESSEDEAG PWILLSDDKG DHPPPPESKI QSFFGLWYRG KAESSEDETS SPAPSKLGGE EEAQPQSPAP DPPCSALHEH LCLGASAAPE A
Source	Baculovirus
Target Names	RILP
Protein Names	Recommended name: Rab-interacting lysosomal protein
Expression Region	1-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.