



Recombinant Human Ankyrin repeat and SOCS box protein 16 (ASB16)

Product Code	CSB-BP850401HU
Abbreviation	ASB16
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96NS5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MARETFPFTS SMLRSLRLQQ EWLEWEDRRR AAAQQCRSRR CPSSPRARLT RPHRSCRDPA VHQALFSGNL QQVQALFQDE EAANMIVETV SNQLAWSAEQ GFWVLTPKTK QTAPLAIATA RGYTDCARHL IRQGAELDAR VGGRAALHEA CARAQFDCVR LLLTFGAKAN VLTEEGTTPL HLCTIPESLQ CAKLLLEAGA TVNLAAGESQ ETPLHVAAR GLEQHVALYL EHGADVGLRT SQGETALNTA CAGAEGPGSC RRRHQAAARRL LEAGADARAA GRKRHTPLHN ACANGCGGLA ELLRYGARA EVPNGAGHTP MDCALQAVQD SPNWEPEVLF AALLDYGAQP VRPEMLKHCA NFPRALEVLL NAYPCVPSCE TWVEAVLPEL WKEHEAFYSS ALCMVNQPRQ LQHLARLAVR ARLGSRQRQ ATRLPLPPLL RDYLLLRVEG CIQ
Source	Baculovirus
Target Names	ASB16
Protein Names	Recommended name: Ankyrin repeat and SOCS box protein 16 Short name= ASB-16
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein is a member of the ankyrin repeat and SOCS box-containing (ASB) family of proteins. They contain ankyrin repeat sequence and a SOCS box domain. The SOCS box serves to couple suppressor of cytokine signalling (SOCS) proteins and their binding partners with the elongin B and C complex, possibly targeting them for degradation.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.