



Recombinant Human Regulation of nuclear pre-mRNA domain-containing protein 1A (RPRD1A)

Product Code	CSB-BP850404HU
Abbreviation	RPRD1A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96P16
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SAFSEAALE KKLSELSNSQ QSVQTLSSLWL IHHRKHSRPI VTVWERELRK AKPNRKLTLFL YLANDVIQNS KRKGPEFTKD FAPVIVEAFK HVSSETDESC KKHLGRVLSI WEERSVYEND VLEQLKQALY GDKKPRKRTY EQIKVDENEN CSSLGSPSEP PQTLDLVRAL QDLENAASGD AAVHQRIASL PVEVQEVSL DKITDKESGE RLSKMVEDAC MLLADYNGRL AAEIDDRKQL TRMLADFLRC QKEALAEKEH KLEEYKRKLA RVSLVRKELR SRIQSLPDLS RLPNVTGSHM HLPFAGDIYS ED
Source	Baculovirus
Target Names	RPRD1A
Protein Names	Recommended name: Regulation of nuclear pre-mRNA domain-containing protein 1A Alternative name(s): Cyclin-dependent kinase inhibitor 2B-related protein p15INK4B-related protein
Expression Region	2-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.