



Recombinant Human Gasdermin-A (GSDMA)

Product Code	CSB-EP850415HU-B
Abbreviation	GSDMA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96QA5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTMFENVTRA LARQLNPRGD LTPLDSLIDF KRFHPFCLVL RKRKSTLFWG ARYVRTDYTL LDVLEPGSSP SDPTDTGNFG FKNMLDTRVE GDVDVPKTVK VKGTAGLSQN STLEVQTL SV APKALETVQE RKLAAHPFL KEMQDQGENL YVMEVTV QEVTLERAGK AEACFSLPFF APLGLQGSIN HKEAVTIPKG CVLAFRVRQL MVKKGKDEWDI PHICNDNMQT FPPGEKSGEE KVILIQASDV GDVHEGFRTL KEEVQRETQQ VEKLSRVGQS SLLSSLSKLL GKKKELQDLE LALEGALDKG HEVTLEALPK DVLLSKEAVG AILYFVGALT ELSEAQQKLL VKSMEKKILP VQLKLVESTM EQNFLDKEG VFPLQPELLS SLGDEELTLT EALVGLSGLE VQRSGPQYMW DPDTLPRLCA LYAGLSLLQQ LTKAS
Source	E.coli
Target Names	GSDMA
Protein Names	Recommended name: Gasdermin-A Alternative name(s): Gasdermin-1
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.