



Recombinant Human Fer3-like protein (FERD3L)

Product Code	CSB-EP850424HU-B
Abbreviation	FERD3L
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96RJ6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAYPESCVD TTVLDFVADL SLASPRRPLL CDFAPGVSLG DPALALREGR PRRMARFEEG DPEEEEECEVD QGDGEEEEEEE ERGRGVSLLG RPKRKRIVITY AQRQAANIRE RKRMFNLNEA FDQLRRKVPT FAYEKRLSRI ETLRLAIVYI SFMTELLESC EKKEG
Source	E.coli
Target Names	FERD3L
Protein Names	Recommended name: Fer3-like protein Alternative name(s): Basic helix-loop-helix protein N-twist Class A basic helix-loop-helix protein 31 Short name= bHLHa31 Nephew of atonal 3 Neuronal twist
Expression Region	1-166
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.