



Recombinant Human Ubiquitin-like protein 7 (UBL7)

Product Code	CSB-BP850429HU
Abbreviation	UBL7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96S82
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSLSDWHLAV KLADQPLTPK SILRLPETEL GEYSLGGYSI SFLKQLIAGK LQESVPDPEL IDLIYCGRKL KDDQTLDFYG IQPGSTVHVL RKSWEPEPDQK PEPVDKVAAM REFRVLHTAL HSSSSYREAV FKMLSNKESL DQIIVATPGL SSDPALGVL QDKDLFSVFA DPNMLDTLVP AHPALVNAIV LVLHSVAGSA PMPGTDSSSR SMPSSSYRDM PGGFLFEGLS DEDDFHPNT RSTPSSSTPS SRPASLGYSG AAGPRPITQS ELATALALAS TPESSSHTPT PGTQGHSSGT SPMSSGVQSG TPITNDLFSQ ALQHALQASG QPSLQSQWQP QLQQLRDMGI QDDELSLRAL QATGGDIQAA LELIFAGGAP
Source	Baculovirus
Target Names	UBL7
Protein Names	Recommended name: Ubiquitin-like protein 7 Alternative name(s): Bone marrow stromal cell ubiquitin-like protein Short name= BMSC-UbP Ubiquitin-like protein SB132
Expression Region	1-380
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.