



Recombinant Arabidopsis thaliana GDSL esterase/lipase 7 (GLIP7)

Product Code	CSB-EP850686DOA
Abbreviation	GLIP7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LFJ9
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	Q SRDHQPLAPA FFVFGDSLVD SGNNNYIPTL ARANYFPYGI DFGFPTGRFC NGRTVVVDYGA TYLGLPLVPP YLSPLSIGQN ALRGVNYASA AAGILDETGR HYGARTTFNG QISQFEITIE LRLRRFFQNP ADLRKYLAKS IIGINIGSND YINNYLMPER YSTSQTYSGE DYADLLIKTL SAQISRLYNL GARKMVLGAS GPLGCIPSQL SMVTGNNTSG CVTKINNMVS MFNSRLKDLA NTLNTTLPGS FFVYQNVFDL FHDMVVNPSR YGLVVSNEAC CGNGRYGGAL TCLPLQQPCL DRNQYVFWDA FHPTETANKI IAHNTFSKSA NYSYPISVYE LAKL
Source	E.coli
Target Names	GLIP7
Protein Names	Recommended name: GDSL esterase/lipase 7 EC= 3.1.1.- Alternative name(s): Extracellular lipase 7
Expression Region	20-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.