



Recombinant Arabidopsis thaliana Dehydration-responsive element-binding protein 2C (DREB2C)

Product Code	CSB-EP850687DOA-B
Abbreviation	DREB2C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LFR2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MPSEIVDRKR KSRGTRDVAE ILRQWREYNE QIEAESCIDG GPKSIRKPP PKGSRKGCMK GKGGPENGIC DYRGVRRRW GKWVAEIREP DGGARLWLGT FSSSYEAALA YDEAAKAIYG QSARLNLPEI TNRSSSTAAT ATVSGSVTAF SDESEVCARE DTNASSGFGQ VKLEDSCSEY VLLDSSQCIK EELKGKEEVR EEHNLA VGFG IGQDSKRETL DAWLMGNGNE QEPLEFGVDE TFDINELGI LNDNNVSGQE TMQYQVDRHP NFSYQTQFPN SNLLGSLNPM EIAQPGVDYG CPYVQPSDME NYGIDLDHRR FNDLDIQDLD FGGDKDVHGS T
Source	E.coli
Target Names	DREB2C
Protein Names	Recommended name: Dehydration-responsive element-binding protein 2C Short name= Protein DREB2C
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.