



Recombinant *Oryza sativa* subsp. japonica Putative ataxin-3 homolog (Os01g0851400, LOC_Os01g63250)

Product Code	CSB-MP850702OFG
Abbreviation	Os01g0851400, LOC_Os01g63250
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LQ36
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MEEAAAASNG GLLYHEVQEG KLCAVHCVNT TLQGPFSEF DLSALAVDLD QRERQVMSEG AAGAATTAAG DFLAEGEGSH NVSLGGDFSI QVLQKALEVW DLQVIPLDSP DVGSCCLFDPE LETAFICHLQ DHWFCIRKVN GEWYNFNSLY PAPEHLSKFY LSAFIDTLKG SGWSIFAVRG NFPKECPMAT EGSNGFGQWL TPDDARRITS SCNQVQPTQ QAGVSLVADQ SEEMSEMDMI AAQQEEADLN AAIAASLMDT GGPFANYAAH EESRSQDAFA IESTSGEMSK DGNLEEQQAN KSETSEPNSD NIESASG SNP KQNTTSLEGK ESIKED
Source	Mammalian cell
Target Names	Os01g0851400
Protein Names	Recommended name: Putative ataxin-3 homolog EC= 3.4.19.12
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.