



Recombinant Mitogen-activated protein kinase pmk-2 (pmk-2)

Product Code	CSB-MP850767CXY
Abbreviation	pmk-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8MXI4
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MGMSATMGDS ASIPGVFFAD FGPAPPEITP EGYHEVELNK TKWVLPQWYN SLKPLGEGAY GVVCTAEYEP TGDRVAIKKF FRPFQSTIHA KRTYRELKLL RTLQHDNVLE MIDVFTDPDP ASSLNNVYFV SVLMGSDLQN IMKIQLTDE QIQLLIYQVL RGLKYIHSAG IHRDLKPSN IAVNERCEVK VFLSFSQLSF LILSFFKILD FGLARAQDAE MTGYVATRWY RAPEIMLNWM HYTQTVDVWS VGCILAEVLS GRPLFPGDDH IDQLTKIMSV VGTPKEEFWS KIQSEARNY IKNRSPIIRQ DFVTLFPMAS PYALELLEMM LILDPDRRIS VSSALRHDYL REYSVPNDEP VAMDTVINSI VTIDPAEERA TTLSDWRELI WNEIRLFQNS ARRLSFVSCT DTEEEPMMKI
Source	Mammalian cell
Target Names	pmk-2
Protein Names	Recommended name: Mitogen-activated protein kinase pmk-2 EC= 2.7.11.24 Alternative name(s): Stress-activated protein kinase pmk-2 p38 MAP kinase 2
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.