



Recombinant Human EP300-interacting inhibitor of differentiation 3 (EID3)

Product Code	CSB-MP850775HU
Abbreviation	EID3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N140
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKMDVSVRAA GCDDLSSGE ADVDPKLEL TADEEKCRSI RRQYRQLMYC VRQNREDIVS SANNSLTEAL EEANVLFDGV SRTREAALDA RFLVMASDLG KEKAKQLNSD MNFFNQLAFC DFLFLVGLN WMEGDPDKLS DCDDSIALSF WKAIEKEATS WMVKAETFHF VFGSFKLERS APKPRLEHQK KVRKMEENGN MPTKLQKLDL SSYPEATEKN VERILGLLQT YFRKYPDTPV SYFEFVIDPN SFSRTVENIF YVSFIVRDGF ARIRLDEDRL PILEPMNVNQ MGEENDSSCH GRKQGVISLT LQEWKNIVAA FEISEAMITY SSY
Source	Mammalian cell
Target Names	EID3
Protein Names	Recommended name: EP300-interacting inhibitor of differentiation 3 Short name= EID-3 Alternative name(s): E1A-like inhibitor of differentiation 3 EID-1-like inhibitor of differentiation 3 Non-structural maintenance of chromosomes el
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.