



# Recombinant Human Kelch-like protein 34 (KLHL34)

<b>Product Code</b>	CSB-EP850783HU
<b>Abbreviation</b>	KLHL34
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N239
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSYFLSYCKA HGGALLTGYQ ALRAEGFLCD VTLETEGSEF PAHRLLACS SDYFRALFKS HTQESRARVI HLHVPSAAGL QRLLDFIYTA WLSLSMDTVE DTLEAASYLQ VTEALGLCGR YLERQLAPEN CCFAANVAAR FGLAHTLDAA ERCIVSHLQE LLARGAGPAG LLELNPTSLR AVLGAPDVAR VPEARLLGLA LAWLRQEPTT ERLAHCTELL ERVRFGLVPA DVLRRVYSGS GLVLPARVKG LIIQALNYHT TPSRQPLMQG EQTSIRSPQT RILLVGGRRR REVVIEEVAA PQRAARGQVA APEPEEEEEEE LEEEEEEEEEW ELTQNVVAFD VYNHRWRSLT QLPTPLLGHS VCTAGNFLFV LGGESPSGSA SSPLADDSRV VTAQVHRYDP RFHAWTEVPA MREARAHFWC GAVGERLLAV GGLGAGGEVL ASVEMYDLRR DRWTAAGALP RALHGHAGAV GDRGVVYISG GKAGRGECCA SSLRDLYVLG PEEQVWSKKA PMGTARFGHH MAVLRGAVFA FLGRYEPFSE IERYDPGADQ WTRLRPLPYD RFCYGLAVVE ETALLGGLK WRDSRQVPTR NVVGYDLDLR RWEDIGCALP WAWSGLRCAV LQLAEGGDDE REGEVGEALD LVLG
<b>Source</b>	E.coli
<b>Target Names</b>	KLHL34
<b>Protein Names</b>	Recommended name: Kelch-like protein 34
<b>Expression Region</b>	1-644
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.