



Recombinant Human Uncharacterized protein C1orf51 (C1orf51)

Product Code	CSB-MP850790HU
Abbreviation	C1orf51
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N365
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSPSSVSSY SSYSLSSSFP TSPVNSDFGF PSDSEREDKG AHGPRPDTVG QRGGSRPSPG PIRCRHRSKV SGNQHTPSHP KQRGSASPMA GSGAKRSRDG ELETSLNTQG CTTEGDLLFA QKCKELQGF PPLTDLLNGL KMGRFERGLS SFQQSVAMDR IQRIVGVLQK PQMGERYLGT LLQVEGMLKT WFPQIAAQKS SLGGGKHQLT KHFPSSHSDS AASSPASPME KMDQTQLGHL ALKPKQPWHL TQWPAMNLTW IHTTPICNPP LSSPGTISFS HGPLGTGTGI GVILFLQHG V QPFTHSAPTT PVPPTTASPV IPGEPMKLSG EGPRCYSLPV TLPDWSYTL SPPSLPTLAR KMTIGHREQQ RSHPPVAADA HLLNL
Source	Mammalian cell
Target Names	CIART
Protein Names	Recommended name: Uncharacterized protein C1orf51
Expression Region	1-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.