



Recombinant Human CAP-Gly domain-containing linker protein 4 (CLIP4)

Product Code	CSB-MP850792HU
Abbreviation	CLIP4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N3C7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTIEDLPDFP LEGNPLFGRY PFIFSASDTP VIFSISAAPM PSDCEFSFFD PNDASCQEIL FDPKTSVSEL FAILRQWVPQ VQQNIDIIGN EILKRGCVN DRDGLTDMTL LHYTCKSGAH GIGDVETAVK FATQLIDLGA DISLRSRWTN MNALHYAAYF DVPELIRVIL KTSKPKDVDA TCSDFNFGTA LHIAAYNLCA GAVKCLLEQG ANPAFRNDKG QIPADVVPDP VDMPLMADA AATAKEIKQM LLDAVPLSCN ISKAMPLNYD HVTGKAMLTSL LGLKLGDRVV IAGQKVGTLR FCGTTEFASG QWAGIELDEP EGKNNGSVGK VQYFKCAPKY GIFAPLSKIS KAKGRRKNIT HTPSTKAAVP LIRSQKIDVA HVTSKVNTGL MTSKKDSASE STLSLPPGEE LKTVTEKDVA LLGSVSSCSS TSSLEHRQSY PKKQNAISSN KKTMSKSPSL SSRASAGLNS SATSTANNSR CEGELRLGER VLVVGQRLGT IRFFGTTNFA PGYWYGIELE KPHGKNDGSV GGVQYFSCSP RYGIFAPPSR VQRVTDSLDT LSEISSNKQN HSYPGFRRSF STTSASSQKE INRRNAFSKS KAALRRSWSS TPTAGGIEGS VKLHEGSQVL LTSSNEMGTV RYVGPTDFAS GIWLGLELRS AKGKNDGSVG DKRYFTCKPN HGVLVPRSRV TYRGINGSKL VDENC
Source	Mammalian cell
Target Names	CLIP4
Protein Names	Recommended name: CAP-Gly domain-containing linker protein 4 Alternative name(s): Restin-like protein 2
Expression Region	1-705
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.