



Recombinant Human Arrestin domain-containing protein 1 (ARRDC1)

Product Code	CSB-BP850807HU
Abbreviation	ARRDC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N5I2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MGRVQLFEIS LSHGRVVYSP GEPLAGTVRV RLGAPLPFRA IRVTCIGSCG VSNKANDTAW VVEEGYFNSS LSLADKGLP AGEHSFPFQF LLPATAPTSF EGPFGKIVHQ VRAAIHTPRF SKDHKCSLVF YILSPLNLNS IPDIEQPNVA SATKKFSYKL VKTGSVVLTA STDLRGYVVG QALQLHADVE NQSGKDTSPV VASLLQKVSY KAKRWIHDVR TIAEVEGAGV KAWRRAQWHE QILVPALPQS ALPGCSLIHI DYYLQVSLKA PEATVTLPVF IGNIAVNHAP VSPRPGLGLP PGAPPLVVPS APPQEEAEAE AAAGGPHFLD PVFLSTKSHS QRQPLLATLS SVPGAPEPCP QDGSPASHPL HPPLCISTGA TVPYFAEGSG GPVPTTSTLI LPPEYSSWGY PYEAPPSYEQ SCGGVEPSLT PES
Source	Baculovirus
Target Names	ARRDC1
Protein Names	Recommended name: Arrestin domain-containing protein 1
Expression Region	1-433
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.