



Recombinant Human Probable allantoicase (ALLC)

| | |
|--------------------------|---|
| Product Code | CSB-YP850815HU |
| Abbreviation | ALLC |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q8N6M5 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MADAPKEGRL TRFLDFTQLM DMASESVGGK ILFATDDFFA PAENLIKSDS PCFKEHEYTE FGKWMGDWET RKRIPGHDW CVLRLGIQGV IRGFDVDVSY FTGDYAPRVS IQAANLEEDK LPEIPERGTR TGAAATPEEF EAI AELKSDD WSYLVPMTEL KPGNPASGHN YFLVNSQQRW THIRLNIFPD GGIARLRVFG TGQKDWTATD PKEPADLVAI AFGGVCVGF S NAKFGHPNNI IGVGGAKSMA DGWETARRLD RPPILENDEN GILLVPGCEW AVFRLAHPGV ITRIEIDTKY FEGNAPDSCK VDG CILTTQE EEAVIRQKWI LPAHKWKPLL PVTKLSPNQS HLFDSL TLEL QDVITHARLT IVPDGGVSRL RLRGFPSSIC LLRPREKPML KFSVSFKANP |
| Source | Yeast |
| Target Names | ALLC |
| Protein Names | Recommended name: Probable allantoicase EC= 3.5.3.4 Alternative name(s): Allantoate amidinohydrolase |
| Expression Region | 1-410 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.