



Recombinant Human APC membrane recruitment protein 2 (FAM123A)

Product Code	CSB-YP850820HU
Abbreviation	FAM123A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N7J2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	METSRSRGGG GAVSERGGAG ASVGVCRRKA EAGAGTGTLA ADMDLHCDCA AETPAAEPPS GKINKAAFKL FKRRKSGGTM PSIFGVKNKG DGKSSGPTGL VRSRTHDGLA EVLVLESGRK EEPGGGGDSG GGGGGRPNPG PPRAAGPGGG SLASSSVAKS HSFFSLLKKN GRSENGKGEP VDASKAGGKQ KRGLRGLFSG MRWHRKDKRA KAEAAEGRAP GGGLILPGSL TASLECVKEE TPRAAREPEE PSQDAPRDPA GEPAGGEEVP APADRAPARS CREAEGLAHP GDTGARGEDA AGHRRRAEPGP GEVRTAEDAS RTGAVPVKTV PLVDSEGGSG RAPAAPDPAS VDPPSDPSAD RICLMFSDVT SLKSFDSL TG CGDIIADQEE EAGPSCDKHV PGPGKPALSK KNPGVVAYQG GGEEMASPDE VDDTYLQEFW DMLSQTEEQG PEPQEGAAGV AAALETKVVP ETPKDTRC VE AAKDASSVKR RRLNRIPIEP HPKEEPKHPE KEQQEGVPNS DEGYWDSTTP GPEEDSSSSG KKAGIPRDSY SGDALYDLYA DPDGSPATLP GGKDNEETSS LSRLKPVSPG TITCPLRTPG SLLKDSKIPI SIKHLTNLPS SHPVVHQQPS RSEMPRTKIP VSKVLVRRVS NRGLAGTTIR ATACHDSAKK L
Source	Yeast
Target Names	AMER2
Protein Names	Recommended name: APC membrane recruitment protein 2 Short name= Amer2 Alternative name(s): Protein FAM123A
Expression Region	1-671
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.