



Recombinant Human FERM domain-containing protein 1 (FRMD1)

Product Code	CSB-MP850827HU
Abbreviation	FRMD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N878
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAVPPRGRGI DPARTNPDTF PPSGARCMEP SPERPACSQQ EPTLGMDAMA SEHRDVLVLL PSREQLRLAV GVKATGRELF QQVCNVASIR DAQFFGLCVV RNNEYIFMDL EQKLSKYFSK DWKKERNEGN EKPRAPFVAF LRVQHYVENG RVIDDHRARH LYYCHLKERV LRSQCAHREE AYFLLAACAL QADLGEHRES AHAGRYFEPH SYFPQWIITK RGIDYILRHM PTLHRERQGL SPKEAMLCFI QEACRLEDVP VHFFRLHKDK KEGRPTVILG LALRGVHIYQ GKKLEIQLDG LPAAQKLVYY TGCTWRSRHL LHLLRASHQL HLRVRPTLQQ LRQREEAEEK QHYRESYISD ELELDLASRS FPGSGVSSQH CPHCLSRHSA DSHGSSYTSG IKANSWLRES REMSVDVPLE VHGLHEKEPS SSPRTSRSH STRGDSQATR QEPCTQVRTR GQSAEAVHQI QEMTAGVSEE QHSHGLDDMQ LHQLALHPAP TSLSHTFHRA LDCRLAGPCE TRATLPSKRS SNCLALDLFG EAPPQEFVV
Source	Mammalian cell
Target Names	FRMD1
Protein Names	Recommended name: FERM domain-containing protein 1
Expression Region	1-549
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.