



Recombinant Human Dynein assembly factor 3, axonemal (DNAAF3)

Product Code	CSB-MP850845HU
Abbreviation	DNAAF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N9W5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> MTPAGSGSG FGSVSWWGLS PALDLQAESP PVDPDSQADT VHSNPELDVL LLGSVDGRHL LRTLSRAKFW PRRRFNFFVL ENNLEAVARH MLIFSLALEE PEKMG LQERS ETFLEVWGNA LLRPPVAAFV RAQADLLAHL VPEPDRLEEQ LPWLSLRALK FRERDALEAV FRFWAGGEKG PQAFPM SRLW DSRLRH YLGS RYDARRGVSD WDLRMK LHDR GAQVIHPQEF RRWRDTG VAF ELRDSSAYHV PNRTLASGRL LSYRGERVAA RGYWGD IATG PFVAFGIEAD DESLLRTSNG QPVKTAGEIT QHNVTELLRD VAAWGRARAT GGDLEEQQHA EGSPEPGTPA APTPESFTVH FLPLDSAQTL HHKSCYNGRF QLLYVACGMV HLLIPELGAC VAPGGNLIVE LARYLVDVRQ EQLQGFNTRV RELAQ AAGFA PQTGARPSET FARFCKSQES ALGNTVPAVE PGTPPLDILA QPLEASNPAL EGLTQPLQGG TPHCEPCQLP SESPGLSEV LAQPQ GALAP PNCESDSK TG V </p>
Source	Mammalian cell
Target Names	DNAAF3
Protein Names	Recommended name: Dynein assembly factor 3, axonemal
Expression Region	1-541
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.