



Recombinant Human DDB1- and CUL4-associated factor 4-like protein 2 (DCAF4L2)

Product Code	CSB-EP850847HU
Abbreviation	DCAF4L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NA75
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MESKRPRLLE EADKQKKTVR VGLNAPSMLR KNQLGFLRFA NYCRIARELR VSCMQRKKVQ IHSWDPSSLA SDRFNRIAN TNTDQLFTVN QVEAGGSKYG IITMRGLTTP ELRVYPHKTLYVPNRKVNMS CWASLNHLDS HLLLCFVGLA DTPSCAVLLP ASLFIGSFPG MRRPGMLCSF QIPDAWSCAW SLSIHAYHSF STGLSQQVLL TNVVTGHQQS FGTSSDVLAQ QFAIMTPLLF NGCRSGEIFG IDLRCGNQGS GWKAICLSHD SAVTSLQILQ DGQFLVSSDM TGTIKLWDLR ATKCVTQYEG HVNNSAYLPV HVNEEEGVVA AVGQDCYTRI WSLRHGHLIT TIPSPYPASE NDIPSVAFSS RLGGFRGAPG LLMVREDLY CFSYG
Source	E.coli
Target Names	DCAF4L2
Protein Names	Recommended name: DDB1- and CUL4-associated factor 4-like protein 2 Alternative name(s): WD repeat-containing protein 21C
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.