



# Recombinant Human NHL repeat-containing protein 2 (NHLRC2)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP850853HU-B   |
| <b>Abbreviation</b>      | NHLRC2   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q8NBF2   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Homo sapiens (Human)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MAAPGGRGRS LSGLLPAQTS LEYALLDAVT QQEKDSLVIYQ YLQKVDGWEEQ<br>DLSVPEFPEG LEWLNTEEPI SVYKDLGKI VVLDFFTYCC INCIHLLPDL<br>HALEHTYSDK DGLLIIGVHS AKFPNEKVLD NIKSAVLRYN ITHPMVNDAD<br>ASLWQELEVS CWPTLVILGP RGNMLFSLIG EGHKDKFLY TSIALKYYKD<br>RGQIRDNKIG IKLYKDSLPP SPLLFPGKVT VDQVTDRLVI ADTGHRILV<br>VWKNQIQYS IGGPNPGRKD GIFSESTFNS PQGVAIMNNI IYVADTENHL<br>IRKIDLEAEK VSTVAGIGIQ GTDKEGGAKG EQQPISSPWD VVFGTSGSEV<br>QRGDILWIAM AGTHQIWALL LDSGKLPKN ELTKGTCLRF AGSGNEENRN<br>NAYPHKAGFA QPSGLSLASE DPWSCLFVAD SESSTVRTVS LKDGAVKHLV<br>GGERDPMNLF AFGDVDGVGI NAKLQHPLGV TWDKRNLLY VADSYNHKIK<br>VVDPKTKNCT TLAGTGDTNN VTSSSFTEST FNEPGGLCIG ENGELLYVAD<br>TNNHQIKVMD LETKMVSVLP IFRSENAVVD GPFLVEKQKT LPKLPKSAPS<br>IRLSPVTACA GQTLQFKLRL DLPSGSKLTE GVSSCWFLTA EGNEWLLQGG<br>IAAGDIENIS SQPTISLQIP DDCLSLEAIV SVSVFLYYCS ADSSACMMKA<br>ILFSQPLQIT DTQQGCIAPV ELRYVF |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | NHLRC2   |
| <b>Protein Names</b>     | Recommended name: NHL repeat-containing protein 2  |
| <b>Expression Region</b> | 1-726  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final   |



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.