



Recombinant Human NHL repeat-containing protein 2 (NHLRC2)

Product Code	CSB-MP850853HU
Abbreviation	NHLRC2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NBF2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAPGGRGRS LSGLLPAQTS LEYALLDAVT QQEKDSLVIYQ YLQKVDGWEEQ DLSVPEFPEG LEWLNTEEPI SVYKDLGKI VVLDFFTYCC INCIHLLPDL HALEHTYSDK DGLLIIGVHS AKFPNEKVLD NIKSAVLRYN ITHPMVNDAD ASLWQELEVS CWPTLVILGP RGNMLFSLIG EGHKDKFLY TSIALKYYKD RGQIRDNKIG IKLYKDSLPP SPLLFPQKVT VDQVTDRLVI ADTGHRILV VWKNQIQYS IGGPNPGRKD GIFSESTFNS PQGVAIMNNI IYVADTENHL IRKIDLEAEK VSTVAGIGIQ GTDKEGGAKG EQQPISSPWD VVFGTSGSEV QRGDILWIAM AGTHQIWALL LDSGKLPKN ELTKGTCLRF AGSGNEENRN NAYPHKAGFA QPSGLSLASE DPWSCLFVAD SESSTVRTVS LKDGAVKHLV GGERDPMNLF AFGDVDGVGI NAKLQHPLGV TWDKRNLLY VADSYNHKIK VVDPKTKNCT TLAGTGDTNN VTSSSFTEST FNEPGGLCIG ENGELLYVAD TNNHQIKVMD LETKMVSVLP IFRSENAVVD GPFLVEKQKT LPKLPKSAPS IRLSPVTACA GQTLQFKLRL DLPSGSKLTE GVSSCWFLTA EGNEWLLQGG IAAGDIENIS SQPTISLQIP DDCLSLEAIV SVSVFLYYCS ADSSACMMKA ILFSQPLQIT DTQQGCIAPV ELRYVF
Source	Mammalian cell
Target Names	NHLRC2
Protein Names	Recommended name: NHL repeat-containing protein 2
Expression Region	1-726
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.