



Recombinant Human Arrestin domain-containing protein 4 (ARRDC4)

Product Code	CSB-MP850861HU
Abbreviation	ARRDC4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NCT1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGGEAGCAAA VGAEGRVKSL GLVFEDERKG CYSSGETVAG HVLLEASEPV ALRALRLEAQ GRATAAWGPS TCPRASASTA ALAVFSEVEY LNVRLSLREP PAGEGIILLQ PGKHEFPFRF QLPSEPLVTS FTGKYGSIQY CVRAVLERPK VPDQSVKREL QVVSHVDVNT PALLTPVLKT QEKMVGCWFF TSGPVLSLAK IERKGYCNGE APIYAEIEN CSSRLIVPKA AIFQTQTYLA SGKTKTIRHM VANVRGNHIA SGSTDTWNGK TLKIPVTPS ILDCCIIRVD YSLAVYIHIP GAKKLMLELP LVIGTIPYNG FGSRNSSIAS QFSMDMSWLT LTLPEQPEAP PNYADVSEE EFSRHIPPYP QPPNCEGEVC CPVFACIQEF RFQPPPLYSE VDPHPSDVEE SQPVSFIL
Source	Mammalian cell
Target Names	ARRDC4
Protein Names	Recommended name: Arrestin domain-containing protein 4
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.