



Recombinant Human Protein FAM13C (FAM13C)

Product Code	CSB-BP850868HU
Abbreviation	FAM13C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NE31
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MFSCFCFSLQ DNSFSSTTVT ECDEDVPSLH EDQTDCCSLR DENNKENYPD AGALVEEHAP PSWEPQQQNV EATVLVDSVL RPSMGNFKSR KPKSIFKAES GRSHGESQET EHVVSSQSEC QVRAGTPAHE SPQNNAFKCQ ETVRLQPRID QRTAISP KDA FETRQDLNEE EAAQVHGVKD PASTQSVL ADGTDSADPS PVHKDGGQNEA DSAPEDLHSV GTSRLLYHIT DGDNPLLSR CSIFSQSQR NLDPEAPSP PSTQQFMMPR SSSRCSCGDG KEPQTITQLT KHIQSLKRKI RKFEKFEQE KKYRPSHGDK TSNPEVLKWM NDLAKGRKQL KELKLLSEE QGSAPKGPPR NLLCEQPTVP RENGKPEAAG PEPSSSGEET PDAALTCLKE RREQLPPQED SKVTKQDKNL IKPLYDRYRI IKQILSTPSL IPTIQEEEDS DEDRPQGSQQ PSLADPASHL PVGDHLTYSN ETEPVRALLP DEKKEVKPPA LSMSNLHEAT MPVLLDHLRE TRADKKRLRK ALREFEEQFF KQTGRSPQKE DRIPMADEYY EYKHIKAKLR LLEVLISKQD VAKTI
Source	Baculovirus
Target Names	FAM13C
Protein Names	Recommended name: Protein FAM13C
Expression Region	1-585
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.