



Recombinant Human Protein FAM83F (FAM83F)

Product Code	CSB-EP850871HU-B
Abbreviation	FAM83F
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEG4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AESQLNCLD EAHVNEKVTE AQAIFYCER RRAALEALLG GGEQAYRERL KEEQLRDFLS SPERQALRAA WSPYEDAVPA ANARGKSKAK AKAPAPAPAE SGESLAYWPD RSDTEVPPLD LGWTDGTFYR GVSRVTLFTH PPKDEKAPHL KQVVRQMIQQ AQKVIADVMD LFTDGDIFQD IVDAACKRRV PVYIILDEAG VKYFLEMCQD LQLTDFRIRN IRVRSVTGVG FYMPMGRIKG TLSSRFLMVD GDKVATGSYR FTWSSSHVDR NLLLLLTGQN VEPFDTEFRE LYAISEEVDL YRQLSLAGRV GLHYSSTVAR KLINPKYALV SGRHPPGEM MRWAARQQRE AGGNPEGQEE GASGGESAWR LESFLKDLVT VEQVLPVPEP IPLGELSQKD GRMVSHMHRD LKPKSREAPS RNGMGAAARG EAAPARRFSS RLFSRRAKRP AAPNGMASSV STETSEVEFL TGKRPNENSS ADISGKTSPS SAKPSNCVIS
Source	E.coli
Target Names	FAM83F
Protein Names	Recommended name: Protein FAM83F
Expression Region	2-500
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.