



Recombinant Human Uncharacterized protein C17orf47 (C17orf47)

Product Code	CSB-MP850873HU
Abbreviation	C17orf47
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEP4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVKTNKPGAK VAVSAQRGSE VTTNTSPQQG HGYVLASSHR SAAVSLNPSH RRSEAAHPTT PHSASDYPRS VSLQSGPGHY AVPTPRGPET GPRTESSRHS SPHLKSQKTQ TLASHASSRQ WKVSPPREEA ARRGESEKSG REVGHHASSI PDAKSTHQLS FQDQKNNLQS QILEDDPPSK VQNPQGVRVP RRILSYPKDE AVQTEPIQRI TTTSEIRSPR SPSLLEHGSS CVSADYQTAQ RRVPEESET GPYGIPIPSKP KALYRNMNLD SLLKLSVLKD SDGVHRVSAR VDPESLHKYS AYPETKPSAK VLVSSQVESN VRTPIRGNSE VGRRVTISPG VQSVEPTHHV TVPSVSEGS KSSMFVTPEP IYKQQTQKPP EITYMSQGPT PRYPELSQKP SIHAELELTP RPLPPRSLPR YGPDSSWWPL LNPEVETPQS QLTPDFEPK CSPSLDLLLS GFKIDSSPFC EDLKFQREKA SLSPSPPKKE FPSWAPLSEV PQTPKHTCKQ PIQRFTAFFL DVSEEMYNRV IWWLKGLCFS LLWAHCGSLG DGRTGEEWHL CIYRAGSFRR
Source	Mammalian cell
Target Names	C17orf47
Protein Names	Recommended name: Uncharacterized protein C17orf47
Expression Region	1-570
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.