



# Recombinant Arabidopsis thaliana Probable protein phosphatase 2C 2 (At1g07160)

<b>Product Code</b>	CSB-BP851401DOA
<b>Abbreviation</b>	At1g07160
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8RX37
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSSSVAVCNS PVFSPSSSLF CNKPLNTSPA HETLTLSLSH LNPPVSSTSP SAASPTSPFC LRLKPPAKL GFGSDSGPGS ILKRKRPTTL DIPVAPV GIA APISNADTPR EESRAVEREG DGYSVYCKRG KREAMEDRFS AITNLQGDPK QAIFGVYDGH GGPTAAEF AA KNLC SNILGE IVGGRNESKI EEAVKRGYLA TDSEFLKEKN VKGGSCCVTA LISDGNLVVA NAGDCRAVLS VGGFAEALTS DHRPSRDDER NRIESSGGYV DTFNSVWRIQ GSLAVSRGIG DAHLKQWIIS EPEINILRIN PQHEFLILAS DGLWDKVSNQ EAVDIARPFC KGTDQKRKPL LACKKLV DLS VSRGSLDDIS VMLIQLCHLF
<b>Source</b>	Baculovirus
<b>Target Names</b>	At1g07160
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 2 Short name= AtPP2C02 EC= 3.1.3.16 Alternative name(s): Protein phosphatase AP2C2
<b>Expression Region</b>	1-380
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.