



# Recombinant Arabidopsis thaliana DNA repair protein recA homolog 2, mitochondrial (At3g10140)

<b>Product Code</b>	CSB-EP851408DOA-B
<b>Abbreviation</b>	At3g10140
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8RY99
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FSYLS QLNGRRSVLA CSGYENRYLS SLVEASDCEL DEVPDDRKVA EKDTALHLAL SQLSGDFDKD SKLSLQRFYR KRRVSVISTG SLNLDLALGV GGLPKGRMVE VYGKEASGKT TLALHIIKEA QKLGGYCAYL DAENAMDPSL AESIGVNTEE LLISRPSSAE KMLNIVDVLV KSGSVDVIVV DSVAAALAPQC ELDAPVGERY RDTQSRIMTQ ALRKIHYSVG YSQT LIVFLN QVRSHVKSNM HFPHAEVTC GGNALPFHAA IRLKMIRTGL IKTANKISGL NVCVQVVKNK LAPGKKKSEL GIHFGHGFYV EREVLELACE HGVILREGTS YFIEGEVIEG KDAAEKYLVE NKEALDTVVA ILRNQLFKM
<b>Source</b>	E.coli
<b>Target Names</b>	At3g10140
<b>Protein Names</b>	Recommended name: DNA repair protein recA homolog 2, mitochondrial Alternative name(s): Recombinase A homolog 2
<b>Expression Region</b>	16-389
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.