



Recombinant *Oryza sativa* subsp. japonica AP2/ERF and B3 domain-containing protein Os01g0693400 (Os01g0693400, LOC_Os01g49830)

Product Code	CSB-BP851411OFG
Abbreviation	Os01g0693400, LOC_Os01g49830
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8RZX9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDSSSCLVDD TNSGGSSTDK LRALAAAAAE TAPLERMGSG ASAVVDAAEP GAEADSGSGG RVCGGGGGGA GGAGGKLPSS KFKGVVPQPN GRWGAQIYER HQRVWLGTFA GEDDAARAYD VAAQRFRGRD AVTNFRPLAE ADPDAAAELR FLATRSKAEV VDMLRKHTYF DELAQSKRTF AASTPSAATT TASLSNGHLS SPRSPFAPAA ARDHLFDKTV TPSDVGKLN LVIPKQHAEK HFPLQLPSAG GESKGVLLNF EDAAGKVWRF RYSYWNSSQS YVLTGWSRF VKEKGLHAGD VVGFYRSAAS AGDDGKLFID CKLVRSTGAA LASPADQPAP SPVKAVRLFG VDLLTAPAPV EQMAGCKRAR DLAATTPPQA AAFKKQCIEL ALV
Source	Baculovirus
Target Names	Os01g0693400
Protein Names	Recommended name: AP2/ERF and B3 domain-containing protein Os01g0693400
Expression Region	1-393
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.