



Recombinant Human Mucin-7 (MUC7)

Product Code	CSB-MP851529HU
Abbreviation	MUC7
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q8TAX7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	RERDHELHRRHHHQSPK SHFELPHYPG LLAHQKPFIR KSYKCLHKRC RPKLPPSPNN PPKFPNPHQP PKHPDKNSSV VNPTLVATTQ IPSVTFPSAS TKITTLPNVT FLPQNATTIS SRENVNTSSS VATLAPVNSP APQDTTAAPP TPSATTPAPP SSSAPPETTA APPTPSATTQ APPSSSAPPE TTAAPPTPPA TTPAPPSSSA PPETTAAPPT PSATTPAPLS SSAPPETTAV PPTPSATTLD PSSASAPPET TAAPPTPSAT TPAPPSSPAP QETTAAPITT PNSSPTTLAP DTSETSAAPT HQTTSVTTQ TTTTKQPTSA PGQNKISRFL LYMKNLLNRI IDDMVEQ
Source	Mammalian cell
Target Names	MUC7
Protein Names	Recommended name: Mucin-7 Short name= MUC-7 Alternative name(s): Apo-MG2 Salivary mucin-7
Expression Region	23-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>