



Recombinant Arabidopsis thaliana Enoyl-CoA hydratase 2, peroxisomal (ECH2)

Product Code	CSB-YP851889DOA
Abbreviation	ECH2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VYI3
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MATSDSEFNS DLLLAHKLPE TRYTYNERDV AIYALGIGAC GQDAVDSDEL KFVYHRNGQD LIQVLPTFAS LFTLGSLTEG LDLPGFKYDP SLLLHGQQYI EIYRPLPSKA SLINKVSLAG LQDKGKAAIL ELETRSYEEG SGELLCMNRT TVFLRGAGGF SNSSQPFSYK NYPSNQGLAV KIPQRQPLTV CEERTQPSQA LLYRLSGDYN PLHSDPEFAK LAGFPRPILH GLCTLGFAIK AIKCVCKGD PTAVKTISGR FLTTVFPGET LITEMWLEGL RVIIYQTKVKE RNKTVLAGYV DIRGLSSSL
Source	Yeast
Target Names	ECH2
Protein Names	Recommended name: Enoyl-CoA hydratase 2, peroxisomal EC= 4.2.1.119
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.